



# 國 盛 生 化 股 份 有 限 公 司

## COUNTRY BEST BIOTECH CO., LTD

台北市忠孝東路二段 112 號 8 樓  
8Fl, No. 112, Sec. 2, Chung-Hsiao East Road, Taipei, Taiwan

TEL: +886-2-23955636-7  
FAX: +886-2-23930556  
E-mail: cbest@cbestbiotech.com  
<https://www.cbestbiotech.com/>

## Product Specification

Product Name: TiterHigh Gold  
**Product Number: CBTH-02**  
Manufactured by: Country Best Biotech Co., Ltd.  
Lot:  
Manufacture Date:  
Expiry Date:  
Storage Temperature: 2 ~ 8 °C  
Emulsion: Oil-in-water

### TEST

### SPECIFICATION

Appearance (Color)	Light yellow to yellow
Appearance (Form)	Liquid
	with particulate
Appearance (Turbidity)	Clear to Slightly Hazy
Forms emulsion with 0.85% NaCl	Pass

Product Sterilized (Gamma Radiation)

#### Disclaimer

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose.